

## **Diastolic Function: What the Sonographer Needs to Know**

---

Pat Bailey, RDCS, FASE  
Technical Director  
Beaumont Health System

---

---

---

---

---

---

---

---

## **Echocardiographic Assessment of Diastolic Function:**

---

Basic Concepts  
Practical Hints and  
Clinical Application

---

---

---

---

---

---

---

---

## **Basic Concepts**

---

- Optimal Diastolic Function occurs when a compliant left ventricle allows filling from low pressure left atrium.

---

---

---

---

---

---

---

---

## Diastolic Dysfunction

- Elevated filling pressures are the main physiologic consequence of diastolic dysfunction.

---

---

---

---


---

---

---

---

## Assessment of Diastolic Function



---

---

---

---


---

---

---

---

## Dystolic Function Exam



---

---

---

---

---

---

---

---

## Left Atrial Volume

- Significant relationship to LA enlargement and diastolic function
    - Measured in apical 4-chamber and two chamber views using area/length method
    - Measure left atrial area
      - Trace volume just prior to MV opening (end systole)
      - Exclude pulmonary veins and appendage
      - Tracing up to the level of the mitral annulus
- (Index the LA Volume to body surface area)  
 LA Volume (cc) / BSA(m<sup>2</sup>)=LA Volume Index

---

---

---

---

---


---

---

---

## Left Atrial Volume

$$\text{LA Volume} = \frac{0.85 \times A_1 \times A_2}{L}$$




---

---

---

---

---

---

---

---

## Mitral Inflow

- Obtain in apical 4-chamber view
  - Technical tips:
    - Perform CW Doppler to assess peak E and A velocities before applying the PW to ensure that maximal velocities are obtained.
    - Using PW Doppler, the cursor must be parallel with the direction of blood flow.

---

---

---

---

---

---

---

---

## Mitral Inflow

- Technical tips cont:
  - A 1-mm to 3-mm sample volume is then placed between the mitral leaflet tips during diastole to record a crisp velocity profile.
  - Optimize gain and wall filter settings to ensure optimal spectral display.
  - Begin with sweep speeds of 25 to 50 mm/s for the evaluation of respiratory variation. If no variation is seen, increase to 100 mm/sec

---

---

---

---

---


---

---

---

## Mitral Inflow

- Align transducer parallel to flow



Adapted from: Appleton, Jensen, Haffe Oh. Doppler Evaluation of Left and Right Ventricular Diastolic Function: Technical Guide

---

---

---

---


---

---

---

---

## Mitral Inflow Variables



---

---

---

---

---


---

---

---



## Tissue Doppler Imaging (TDI)




---

---

---

---

---

---

---

---

## Tissue Doppler Imaging (TDI)

- Measurements:
  - Mitral annular e' velocity (m/sec)
  - Mitral annular a' velocity (m/sec)
    - Limitation for septal annulus:
      - Regional wall motion abnormalities
      - Prosthetic valves
      - Severe mitral valve disease
      - Constrictive pericarditis

---

---

---

---


---

---

---

---

## Pulmonary Venous Flow



Adapted from: The Echo Manual - Second Edition

---

---

---

---

---

---

---

---

## Pulmonary Veins

- Technical tips:
  - Performed in apical 4-chamber view
  - Color flow imaging is useful for the proper location of the sample volume in the RUPV.
  - Angle transducer superiorly with aortic valve in view.
  - Use a 2-mm to 3-mm sample volume, placed >5 cm into the pulmonary vein.
  - Optimize scale/baseline
  - Increase sample volume size if Doppler signal is weak (4-5mm)

---

---

---

---

---

---

---

---

---

---

## Pulmonary Veins

- Measurements:
  - Systolic flow velocity [S<sub>1</sub>,S<sub>2</sub>] (m/sec)
  - Diastolic flow velocity (m/sec)
  - Atrial reversal duration (msec)
    - Limitations:
      - Wall noise (ensure the atrial reversal flow signal intensity matches that of the systolic and diastolic flow)

---

---

---

---

---

---


---

---

---

---

## Pulmonary Veins




---

---

---

---

---

---


---

---

---

---

## Mitral and P-Vein "A" duration comparisons



### AR Duration

- Optimize atrial reversal Doppler signal
- Requires sharp distinct signals
- Measure close to baseline
- Difference in durations must be > 30 msec

---

---

---

---


---

---

---

---

## Comparison of Durations




---

---

---

---

---

---

---

---

## Comparison of Durations

- MV A duration > P vein AR duration in hearts with normal filling pressure
- As LVEDP increases, the P vein AR duration lengthens
- Difference in durations must be > 30 msec to be significant
- Longer P vein AR duration (> 30 msec) indicates increased LVEDP

---

---

---

---

---

---

---

---

## Valsalva Maneuver

- The Valsalva maneuver is used to differentiate the grades (II, III) of diastolic dysfunction.
  - **Performance of valsalva maneuver:**
    - Ask patient to suspend their breath ant the end of a normal inspiration and strain down without breathing. By increasing the intra-thoracic pressures and thereby reducing the venous return to the heart, a decreased pre-load (left atrial pressure) potentially unveils any underlying grade of diastolic dysfunction.

---

---

---

---


---

---

---

---

## Mitral inflow with Valsalva




---

---

---

---

---

---

---

---

## Valsalva Maneuver

- An adequate Valsalva maneuver is defined as a 10% reduction in the maximal E velocity from baseline.
- A reduction of the E/A ratio by 0.5 or more demonstrates increased filling pressure (ex. 1.2 to 0.7)

---

---

---

---

---

---

---

---

## Valsalva Maneuver

- Technical tips:
  - Obtain apical 4-chamber view
  - Align sample volume at leaflet tips
  - Begin Valsalva maneuver, while maintaining sample volume alignment
  - After 12 seconds, turn to live Doppler to capture the signal.

---

---

---

---

---

---

---

---

## Valsalva Maneuver

- Measurements:
  - MV "E" velocity (m/sec)
  - MV "A" velocity (m/sec)
  - E at A velocity (m/sec)
  - E/A ratio
  - Deceleration time (msec)

---

---

---

---

---

---

---

---

## Valsalva Maneuver

- Limitations:
  - Not everyone is able to perform this maneuver adequately

---

---

---

---


---

---

---

---

## Practical Approach to Grade Diastolic Dysfunction




---

---

---

---

---

---

---

---

---

---

## Grades of Diastolic Dysfunction

- Grade I - Abnormal relaxation
- Grade II - Pseudonormal
- Grade IIIa / IIIb - Restrictive

---

---

---

---

---

---


---

---

---

---

## Abnormal LV Filling Patterns



Adapted from: The Echo Manual - Second Edition

---

---

---

---

---

---

---


---

---

---

## Grade I Abnormal Relaxation

- Cont:
- LA contraction (A wave) contributes up to 35% of filling
  - E/A ratio  $\leq 0.8$  (If E at A is  $< 0.2$  m/sec)
  - DT  $> 240$




---

---

---

---

---


---

---

---

## Grade II Pseudo-normal

- Impaired myocardial relaxation with mild to moderate elevation of LV filling pressures.
- Symptoms of heart failure at rest or with minimal exertion (NYHA II-III)
  - E/A 0.8-1.5
  - DT 160-200 ms
  - Av. E/e' 9-12
  - Ar-A  $\geq 30$  ms
  - Val E/A  $\geq 0.5$




---

---

---

---

---

---

---

---

## Grade II Pseudo-normal

- Pseudonormal cont:
  - Pulmonary vein S/D  $\leq 0.5$
  - LV dysfunction (not always)
  - LA enlargement

---

---

---

---

---

---

---

---

### Grade IIIa - IIIb

- LV is less compliant (Incr. filling pressure, Incr. E)
- Rapid equalization of pressure (decr. DT)
- Late diastolic LA-LV pressure gradient is limited (decr. A)
- Symptoms of heart failure at rest or with minimal exertion (NYHA III-IV)
- Reversible mitral inflow pattern upon preload reduction with Valsalva, nitroglycerin or diuretic administration [Grade IIIa]

---

---

---

---

---

---

---

---

### Grade IIIa - IIIb

- Cont:
- Eventually irreversibly restrictive (no response to aggressive preload reduction) [Grade IIIb]



---

---

---

---

---

---

---

---

### THANK-YOU



---

---

---

---

---

---

---

---